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Roma, 19 luglio 1999

Il Ministro: BINDI

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HIGH PERFORMANCE LIQUID CHROMATOGRAPHIC METHOD FOR THE DETERMINATION OF OCHRATOXIN A IN ROASTED COFFEE

1. SCOPE

This document specifies a method for the determination of ochratoxin A in roasted coffee at concentrations up to 10 µg/kg using immunoaffinity column clean up and high performance liquid chromatography.

2. NORMATIVE REFERENCES

This document incorporates by dated or undated reference, provisions from other publications. These normative references are cited at the appropriate places in the text and the publications are listed hereafter. For dated references, subsequent amendments to, or revisions of any of these publications apply to this document only when incorporated in it by amendment or revision. For undated references the latest edition of the publication referred to applies.

EN ISO 3696 Water of analytical laboratory use - Specification and test methods.

3. PRINCIPLE

Ochratoxin A is extracted from roasted coffee by blending with methanol and sodium bicarbonate. The sample is cleaned up by passing first through a phenyl silane column and then through an immunoaffinity column. Ochratoxin A is separated and identified by reversed phase HPLC and quantified by fluorescence.

WARNING

Ochratoxin A is a potent nephrotoxin and liver toxin and has been reported to have immune-suppressant properties.

Gloves and safety glasses should be worn at all times, and all standard and sample preparation stages should be carried out in a fume cupboard.

Toluene is highly flammable and harmful. Standard preparation involving this solvent must be performed in a fume cupboard. Operations outside the fume cupboard, such as measurement of calibrants by UV spectrophotometry, must be performed with the calibrants in closed containers.

4. REAGENTS

During analysis, unless otherwise stated, use only reagents of recognised analytical grade and only distilled water or water of grade 1 according to EN ISO 3696.

4.1 Ochratoxin A calibrant solution (a 10 µg/mL ochratoxin A calibrant in toluene+acetic acid, (99+1, v+v) is supplied). The concentration of this calibrant may be checked using a UV spectrophotometer ($\epsilon = 5440$). The concentration must be taken as exactly 10 µg/mL for all calculations.

4.2 Methanol (HPLC grade)

4.3 Sodium bicarbonate

4.4 3% sodium bicarbonate solution (w+v)

Dissolve 30 g sodium bicarbonate in 1000 mL water in a 1 litre volumetric flask

4.5 Extraction solvent

Mix 1 part per volume methanol (4.2) with 1 part per volume of 3% sodium bicarbonate solution (4.4)

4.6 Phenyl silane column wash reagent 1:

Mix 25 parts per volume methanol (4.2) with 75 parts per volume 3% sodium bicarbonate solution (4.4)

4.7 Phenyl silane column wash reagent 2:

1% sodium bicarbonate solution

Dissolve 1 g sodium carbonate in 100 mL water in a 100 mL volumetric flask

4.8 Phenyl silane column elution reagent

Mix 7 parts per volume methanol (4.2) with 93 parts per volume water

4.9 Phosphate buffered saline

Dissolve 8 g sodium chloride, 1.16 g disodium hydrogen orthophosphate, 0.2 g potassium dihydrogen phosphate and 0.2 g potassium chloride in 1 litre distilled water.

Adjust pH to 7.4 with alkali

4.10 Acetonitrile (HPLC grade)**4.11 Glacial acetic acid****4.12 HPLC mobile phase**

Mix 102 parts per volume water with 96 parts per volume acetonitrile (4.10) and 2 parts per volume acetic acid (4.11) and de-gas (eg with helium)

4.13 Toluene (AR grade)**4.14 Toluene acetic acid mixture**

Mix 99 parts per volume toluene (4.14) with 1 part per volume acetic acid (4.15)

4.15 Standard dilution solution (Methanol + water, 3+7, v+v)

Mix 3 parts per volume methanol (4.2) with 7 parts per volume water

4.16 Silanising reagent, (e.g. Surfasil - Pierce 4289 is suitable)

5. APPARATUS

Normal laboratory glassware and apparatus, and in particular, the following:

5.1 Silanised glass vials are needed to ensure the stability of ochratoxin A in aqueous solvents. Prepare vials by filling them with silanising reagent (4.10) and leave this reagent in the vial for 1 minute. Rinse the vial twice with a low polarity solvent (eg chloroform) then with methanol. Finally wash the vials twice with distilled water.

5.2 Analytical balance accurate to 2 decimal places

5.3 Wrist action shaker; or similar

5.4 Centrifuge capable of speeds up to 2400 rpm and centrifuging at 4 °C

5.5 Centrifuge tubes, 50 mL capacity

5.6 Displacement pipettes (5 mL, 1 mL, 200 µl) with appropriate pipette tips

5.7 Vacuum manifold to accommodate phenyl silane columns and immunoaffinity columns

5.8 Reservoirs and attachments to fit to immunoaffinity columns

5.9 Vacuum pump - capable of pulling a vacuum of 10 mbar and pumping 18 litres per minute

5.10 Filter papers (eg. Whatman number 4)

5.11 Phenyl silane solid phase extraction columns, 500 mg (e.g. JT Baker 7095-03 is suitable

5.12 HPLC apparatus comprising the following:

5.12.1 Valve injection system with 100 μ L injection loop

5.12.2 Mobile phase pump (isocratic) capable of pumping 1 mL/minute pulse free

5.12.3 HPLC column - C18 reverse phase, ODS 1 or 2, 5 μ m, plus a suitable corresponding reversed-phase guard column

5.12.4 Fluorescence detector with flow cell

Emission wavelength 420 nm

Excitation wavelength 390 nm

5.12.5 Data collection system

5.13 Spectrophotometer (optional for checking concentration of supplied calibrant)

5.14 Immunoaffinity columns specific for ochratoxin A clean-up. The column should have a total capacity of not less than 2700 ng ochratoxin A and should give a recovery of not less than 85% when applied in a calibrant solution of ochratoxin A in methanol/PBS containing 100 ng of ochratoxin A.

5.15 Disposable syringe filters 0.2 μ m pore size, 25 mm polysulfone membrane

5.16 Steam bath

6. PROCEDURE

6.1 Extraction

Weigh, to the nearest 0.1 g, a 15 g test portion of roasted coffee sample into a 500 mL conical flask. Add 150 mL extraction solvent (methanol + 3% sodium bicarbonate solution, 50+50, v+v, 4.5). Stopper the flask and shake gently for 30 minutes. Filter the extract through filter paper (5.10). Remove an approximately 50 mL portion of coffee filtrate into a centrifuge tube (5.5) and centrifuge for 15 minutes, 4 °C, 2400 rpm.

6.2 Phenyl silane column clean-up

Attach a phenyl silane column (5.11) to the vacuum manifold (5.7) but do not apply vacuum to the column. Wash the column with 15 mL methanol (4.2) followed by 5 mL 3% sodium bicarbonate solution (4.4). Discard the washings. Pipette 10 mL centrifuged coffee extract (6.1) into a glass beaker and add 10 mL 3% sodium bicarbonate solution (4.4). Pass this diluted extract through the phenyl silane column at a maximum speed of 5 mL/min. Wash the phenyl silane column with 10 mL wash reagent 1 (methanol + 3% sodium bicarbonate solution (25+75, v+v) (4.6)); followed by 5 mL 1% sodium bicarbonate solution (wash reagent 2) (4.7). Remove the phenyl silane column from the vacuum manifold and place over a suitable collection vessel. Elute ochratoxin A with 10 mL phenyl silane column elution reagent (methanol + water (7+93, v+v) (4.8)). The phenyl silane column should be dried with at least three 10 mL volumes of air using a syringe before eluting the ochratoxin.

6.3 Immunoaffinity column clean-up

Dilute the 10 mL of phenyl silane column eluate with 30 mL phosphate buffered saline (PBS) (4.9). The use of a fine pre-filter before the extract is applied to the column is not necessary. Connect the immunoaffinity column (5.14) to the vacuum manifold (5.7), and attach the reservoir (5.8) to the immunoaffinity column. Add the diluted sample extract to the reservoir and pass through the immunoaffinity column. Flow rate should not exceed 5 mL/min. The immunoaffinity column must not be allowed to run dry. Wash the immunoaffinity column with 10 mL water. Remove the immunoaffinity column from the vacuum manifold and place over a silanised vial (5.1). Elute ochratoxin A into the silanised vial with four 1 mL portions of methanol (4.2).

Evaporate immunoaffinity column eluate to dryness over a steam bath (5.16), under nitrogen. Redissolve in 1 mL of HPLC mobile phase (4.12) which has been filtered through 0.2 μ m filter (5.15). Transfer to HPLC vial.

7. SPIKING PROCEDURE

7.1 Preparation of 200 ng/mL spiking solution

Pipette 1 mL of the 10 μ g/mL ochratoxin A calibrant solution (4.1) into a glass vial and dilute to 5 mL with 4 mL with toluene + acetic acid (99+1, v+v) (4.14). This gives a solution of 2 μ g/mL ochratoxin A.

Pipette 1 mL of the 2 μ g/mL ochratoxin A solution into a silanised glass vial. Evaporate the solvent under a stream of nitrogen. Redissolve in 10 mL methanol (4.2). This gives a solution of 200 ng/mL ochratoxin A.

7.2 Weigh, to the nearest 0.1 g, 15 g blank roasted coffee sample (supplied) into a 500 mL conical flask. Pipette 375 μ L of 200 ng/mL ochratoxin A spiking solution (7.1) onto the blank coffee. Proceed immediately as in section 6.1. Note that the spiked coffee should NOT be left overnight and should only be left for a maximum of 3 minutes before analysis.

8. HPLC DETERMINATION

8.1 Calibration graph

Prepare a calibration graph at the beginning of the analysis and whenever the chromatographic conditions change.

8.1.1 Preparation of calibrants:

Pipette 1 mL of the 10 µg/mL ochratoxin A calibrant solution provided (4.1) into a glass vial and dilute to 5 mL with 4 mL toluene/acetic acid (99+1, v+v) (4.14). This gives 2 µg/mL ochratoxin A solution.

Pipette 1 mL of the 2 µg/mL ochratoxin A solution into a silanised glass vial. Evaporate the solvent under a stream of nitrogen. Redissolve in 10 mL methanol (4.2) which has been filtered through a 0.2 µm filter. This gives a 200 ng/mL solution.

Dilute this with calibrant dilution solution (methanol + water (3+7, v+v) 4.15), to prepare HPLC calibrants.

Prepare five HPLC calibrants in separate 10 mL volumetric flasks according to the table below: Make up each standard to volume (10 mL) with filtered methanol + water (3+7, v+v).

	Cal 1	Cal 2	Cal 3	Cal 4	Cal 5
µL filtered methanol + water (3+7, v+v)	9987.5	9975	9950	9875	9750
µL 200 ng/mL OTA	12.5	25	50	125	250
conc OTA ng/mL	0.25	0.5	1.0	2.5	5.0

Cal :- calibrant ; conc :- concentration

8.2 Operating conditions

When the column specified in 5.12.3 and the mobile phase specified in 4.13 were used, the following settings were found to be appropriate:-

Flow rate: 1 mL/min

Column oven temperature: 35 ± 1°C

Fluorescence detection: Emission wavelength $\lambda = 460$ nm
 Excitation wavelength $\lambda = 333$ nm
 Injection volume: 100 μ L

9. EVALUATION

- 9.1 Determine from the calibration graph, the masses in pg of the ochratoxin A in the aliquot of test solution injected onto the HPLC column.
- 9.2 Calculate the mass fraction, W_{OTA} , of ochratoxin A in micrograms per kilogram using the equation:

$$W_{OTA} = \frac{M_A \times V_4}{V_3} \times \frac{V_1}{V_2} \times \frac{1000}{M_s} \times 10^{-6}$$

Where: M_A = mass of ochratoxin A in aliquot of test solution injected on column in pg (9.1)

V_4 = volume of aliquot of test solution injected on column, in mL (8.2)

V_3 = volume of test solution (1 mL, 6.3)

V_2 = volume of sample filtrate taken for clean-up, in mL (10 mL, 6.2)

V_1 = volume of extraction solvent, in mL (150 mL 6.1)

M_s = mass of sample extracted, in g (15 g 6.1)

10^{-6} = factor to convert pg to μ g

Report the results (on the form supplied) in μ g/kg to one decimal place or more in order to avoid a series of identical results.

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SUMMARY

HPLC determination of benzoic and sorbic acids in derived fruit products

A solid phase extraction technique on a disposable cartridge packed with octadecylbonded silica is proposed as a clean-up step for the determination of benzoic and sorbic acids in derived fruit products (jam samples). After elution from the cartridge the analytes are separated by HPLC on a C18 reversed-phase column and determined using UV absorbance measurements at 234 nm.

Linearity of the response, sensitivity, repeatability of the analytical procedure, recovery studies from spiked jam samples have been carried out. Detection limits are 5 and 3 ppm for benzoic and sorbic acid respectively in apricot and raspberry jam samples.

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9. LEE H.S., ROUSEFF R.L. and FISHER J.F.—Determination of food preservatives in orange juice by reversed phase liquid chromatography.
J. Food Sci., 1986, 51, 568—570.
10. BENNETT M.C. and PETRUS D.R.—Quantitative determination of sorbic acid and sodium benzoate in citrus fruit.
J. Food Sci., 1977, 42, 1220—1221.
11. PUTTEMANS M.L., DRYON L. and MASSART D.L.— Extraction of organic acids by ion-pair formation with tri-n-octylamine. Part V. Simultaneous determination of synthetic dyes, benzoic acid, sorbic acid and saccharin in soft drinks and lemonade syrups.
J. Assoc. Off. Anal. Chem., 1984, 67, 880—885.
12. MURRAY J. — Absorption of benzoic and sorbic acids by herring and prawns in marinades.
International Journal of Food Science and Technology, 1988, 23, 171—175.
13. Milch, Milchpulver, Joghurt und andere gesäuerte Milch, Bestimmung des Gehalts an Benzoe und Sorbinsäure.
Vorläufiger Internationaler IMV Standard 139 (1987).
Milchwissenschaft, 1988, 43, 366—368.
14. NAGAYAMA T., NISHIJIMA M., YASDA K., SAITO K., KAMIMURA H., IBE A., USHIYAMA H., NAGAYAMA M. and NAOI Y.—Benzoic acid in fruits and fruit products.
J. Food Hyg. Soc. Jpn. 1983, 24, 416—422.

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SUMMARY Method of analysis of cyclamate in foods

Sodium and calcium cyclamate are new sweeteners authorized in the E.U. The authors propose a chromatographic method to determine the quantity of this sweetener in food by ion chromatography and conductivity detector. The method was tested on various food products and the repeatability, limit of determination and recovery are specified.

Key-words : *cyclamate, sweetener, ion chromatography, quantity determination.*

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Determination of annatto in high-fat dairy products, margarine and hard candy by solvent extraction followed by high-performance liquid chromatography

FRANK E. LANCASTER and JAMES F. LAWRENCE

Food Research Division, Bureau of Chemical Safety, Food Directorate,
Health Protection Branch, Ottawa, Ontario K1A 0L2, Canada

(Received 16 February 1994; revised 17 March 1994; accepted 11 April 1994)

Utilizing solvents such as ethanolic aqueous ammonia, petroleum ether, hexane and chloroform, annatto components α - and β -norbixin and α - and β -bixin were extracted from cheese, butter, margarine and hard candy. After transferring the extract into a solution of aqueous acetic acid in methanol, bixin and norbixin were determined quantitatively using high-performance liquid chromatography (HPLC) and an absorbance detector set at 500 nm. Recovery of norbixin from spiked cheese samples averaged 92.6% over a range of 1 to 110 $\mu\text{g/g}$. Commercial cheese samples were found to contain 1.1–68.8 $\mu\text{g/g}$ total norbixin, and two samples also contained 5.1–5.6 $\mu\text{g/g}$ total bixin. Samples of uncoloured butter were spiked with bixin and recovery averaged 93.2% over a range of 0.1 to 445 $\mu\text{g/g}$. Levels of 0.2 $\mu\text{g/g}$ total bixin and 0.91 $\mu\text{g/g}$ total norbixin were found in one commercial butter sample; the others contained trace levels of both compounds. Hard candies were prepared in the laboratory and recovery studies conducted. Recovery of norbixin averaged 88%.

Keywords: annatto, bixin, norbixin, food colour, cheese, butter, margarine, hard candy, HPLC

Introduction

Increased public awareness of food colours has produced a stronger demand to substitute natural colours, such as annatto, for synthetics. Natural colours are sometimes perceived, often without valid proof, to be less harmful and more acceptable for use in foods. Annatto, cochineal and β -carotene are three of the most popular natural colours. While the latter has been widely researched, it appears that methodology, particularly quantitative, for the former two is less extensive.

Since the last century, annatto has been a popular colouring agent for cheese and butter. The principal colouring constituents in annatto are carotenoids. Bixin (figure 1) occurs in oil-soluble annatto in *cis* and *trans* forms, known as α - and β -bixin respectively, because of geometrical isomerism around the 5,6-carbon atoms. Water-soluble annatto contains α - and β -norbixin, similar in structure to the bixins, except that the terminal methyl ester group has been removed by saponification, and both carboxyl groups are in their sodium or potassium salt form (Smith *et al.* 1983). Three decades ago, the composition of oil-soluble annatto colours was described by McKeown (1961, 1963, 1965) and McKeown and Mark (1962). Since that time, other researchers have focused on the properties of both bixin and water-soluble norbixin (Reith and Gielen 1971), or the extraction and chemistry of annatto (Preston and Ricard 1980). For many years, annatto compounds have been studied using thin-layer chromatography (TLC) (Kanada *et al.* 1981, Tirimanna 1981, Miyamoto *et al.* 1991), and in 1984 Amakawa *et al.* published what appears to be the first method incorporating high-performance liquid chromatography (HPLC).


Figure 1. Structures of bixin.

In 1991, it was decided to conduct a study of annatto in order to update and modernize our methodology, to review and assess published methods, and if necessary develop modern methods for the determination of annatto in the food group in which it is used most extensively, i.e. high fat dairy products. Water-soluble annatto also appears in confectionery products, snack foods and ice cream. Therefore, it was decided to test first the effects of sugar and heat on water-soluble norbixin before examining more complex food items. Candy-making has proved to be a useful technique (Lancaster and Lawrence 1986, 1992) when examining a food colour for stability because both heat and sugar are utilized. Additional information, for example the effect of acidic ingredients such as citric acid, is also provided. It was also expected that some of the analytical techniques developed during a candy study could be utilized when analysing more complex food products.

Also, commercial samples of butter and cheese were examined, and recovery and degradation studies conducted. The methods described here resulted from that study.

Experimental

Materials

Concentrated (90%) bixin crystals (Warner-Jenkinson Canada) and norbixin powder (containing 28.6% and 13.2% norbixin, and additional components extracted by alkaline media) (Warner-Jenkinson Canada and Pointing Canada, respectively) were used for recovery studies, as spiking standards and during candy preparation. Purity was determined using HPLC and spectrophotometry. Stock solutions of norbixin (150 mg norbixin powder/100 ml 0.2% sodium hydroxide in water) and bixin (50 mg bixin crystals/100 ml 0.08% sodium hydroxide in methanol) were prepared no more than 24 h before use.

All solvents were of HPLC or analytical-reagent (AR) grade. Petroleum ether (30°-60°C) (Caledon) and chloroform (Caledon) were distilled in glass. Celite 545[®] (Fisher) was combined with all samples before extraction.

Preparation of hard candy

A 100 g portion of granulated sugar, 50 g corn syrup, 1 g D-glucose, 1 g citric acid and 50 ml distilled water were combined in a tared 600 ml beaker, and mixed thoroughly. The beaker was transferred to a hot plate (medium heat), a candy thermometer was suspended in the sugar mixture and the solution was heated, and stirred occasionally. When the temperature of the candy reached 135°C, 5 ml water-soluble annatto containing a known amount of norbixin was added. For example, 1.43 g norbixin powder, which was 28.6% norbixin, was combined with 0.2 g NaOH in a 100 ml volumetric flask, diluted to volume with distilled water, and shaken occasionally until the powder dissolved. A 5 ml aliquot of this solution was combined with the hot candy and mixed thoroughly. When the temperature of the candy reached 152°C, the heat source was removed.

The hot candy was poured into a tared aluminium tray, allowed to cool for 15 min, then transferred to a freezer. The beaker containing the residue was allowed to cool, then weighed after it reached room temperature. The difference between this weight and the original weight of the empty beaker was added to the weight of the candy in the tray. After the candy hardened, the tray was removed from the freezer, weighed, and total yield determined (ca 152 g). The candy was transferred to a plastic bag, sealed, and stored in the freezer until required.

Extraction of hard candy

A 20.0 g portion of uncrushed candy was transferred quantitatively into a tared 230 ml glass centrifuge tube, 50.0 ml distilled water was added, and the tube was capped then shaken mechanically for ca 7 min or until the candy dissolved. If spiking was required, a known amount of norbixin was combined with the candy solution. In a fume hood, 50.0 ml 0.5% acetic acid in chloroform was added, and the tube was capped then shaken for 1 min. The cap was loosened then tightened again after 30 s in order to release pressure produced by the chloroform. The solution was centrifuged at 1800 rpm for 10 min.

The chloroform layer was removed using a 25 or 50 ml pipette and bulb, filtered by draining through a chloroform-washed glass wool pledget in the stem of a powder funnel, and collected in a 100 ml volumetric flask. The aqueous layer was retained for further extraction, as described below. The funnel and wool were rinsed with chloroform in order to remove remaining colour, and the wash collected in the same flask. A 0.25 ml portion of acetic acid was added to the flask, the contents diluted to volume with chloroform then mixed. An appropriate volume of extract (e.g. 2.0 ml) was transferred into a 10 ml volumetric flask and diluted to volume with 0.5% acetic acid in chloroform and retained for spectrophotometry or liquid chromatography.

The aqueous layer required two additional extractions. For each, 25.0 ml 0.5% acetic acid in chloroform was added to the candy solution in the centrifuge tube, and the contents shaken and centrifuged, as described above. The chloroform layer was removed and filtered (as above), the extract and wash collected in a 50.0 ml volumetric flask, diluted to volume with 0.5% acetic acid in chloroform and mixed. The concentration of norbixin in each of the three extracts was determined using spectrophotometric curves and absorptivity values of 2620 at 473 nm and 2290 at 502 nm (Smith *et al.* 1983), or the HPLC systems described below.

Extraction of cheese, butter and margarine

Bulk cheese was grated finely before analysis, but process cheese food and butter did not require any preparation. A 20.0 g sample was transferred quantitatively into a 250 ml glass centrifuge tube (tube A) containing 3 g Celite 545 (see also figure 2). In a fume hood, 25 ml of ethanol-water-concentrated ammonium hydroxide (29.7% NH_3), 100:35:15, was combined with the sample. The contents of the tube were mixed well using a glass stirring rod before 50 ml of petroleum ether was added. The tube was capped then shaken for 1 min.

After centrifuging at 1700 rpm for 10 min, the ether layer (top) was removed using a 25 ml pipette (with bulb attached) followed by a Pasteur pipette, and discarded. A second 50 ml portion of petroleum ether was added, the extraction was repeated, and the tube centrifuged before the ether layer was removed and discarded. The ethanolic aqueous ammonia layer was decanted into a second 250 ml centrifuge tube (tube B) and retained.

A second 25 ml portion of ethanolic ammonium hydroxide was added to the Celite/cheese residue in tube A, the contents mixed well, centrifuged, then the liquid was decanted into tube B. The extraction was repeated a third time or until no yellow colour remained on the Celite. A 50 ml portion of 11% aqueous acetic acid was combined with the solution in tube B, mixed well, then 45 ml chloroform-glacial acetic acid (98.5:1.5) were added. Often, proteinaceous residue precipitated. The tube was capped then shaken well for 1 min. During this step the cap was loosened once or twice to relieve pressure, then the tube was centrifuged.

The upper layer (aqueous) was removed and discarded; all annatto transferred to the chloroform. The chloroform layer was filtered through a small pledget of chloroform-washed glass wool in the stem of a powder funnel, and collected in a third 250 ml centrifuge tube (tube C). The glass wool was washed with a minimal amount of chloroform until all colour was removed. A second portion (25 ml) of chloroform-acetic acid was added to the residue in tube B, mixed well and the tube was shaken and centrifuged. The extract was filtered through the glass wool pledget, and collected in tube C. The process was repeated until all colour had been extracted. In most instances, three extractions were sufficient.

A 50 ml portion of 1% aqueous acetic acid was added to the chloroform solution, and the tube was capped and shaken gently for 1 min then centrifuged. The aqueous (top) layer was removed using pipettes, and the chloroform layer was filtered through a second chloroform-washed glass wool pledget in the stem of a powder funnel, collected in a volumetric flask of appropriate size (100–200 ml) and diluted to volume with chloroform. When an emulsion was generated and remained in the funnel, it was washed with an appropriate amount of methanol (25–50 ml), and the wash collected in a separate volumetric flask (25 or 50 ml), diluted to volume with methanol, and retained for HPLC analysis.

A 5.0–10.0 ml aliquot (volume recorded) of the chloroform extract was transferred to a 25 ml Erlenmeyer flask, evaporated to dryness with a stream of nitrogen, then 5.0 ml methanol was added and the colour redissolved. The contents of this solution, the spiking standard, and the methanol wash, were analysed using the HPLC system described below.

Liquid chromatography

The system consisted of a model 110B pump from Beckman, an IBM PS/2 computer system loaded with Beckman GOLD software and connected to a Beckman


Figure 2. Procedure for the extraction of annatto from cheese, butter and margarine.

406 analog interface, a variable wavelength detector (Spectroflow 783; wavelength, 500 nm; detector range, 0.05–0.8 au/s) from Kratos, a Spectra-Physics SP4270 integrator, and a Rheodyne 7725 injection port with 20 μ l loop connected to a Supelco LC-18 (5 μ) column (250 mm \times 4.6 mm i.d.). The mobile phase was methanol–2% aqueous acetic acid (9:1), and flow rate was 1.0 ml/min.

A known amount (5–10 ml) of extract was transferred quantitatively to a 25 ml Erlenmeyer flask, and the solvent was evaporated to dryness with a stream of nitrogen. For most extracts, a 2.0 ml portion of methanol–2% aqueous acetic acid (9:1) was added to the flask, and the colour was redissolved. When examining concentrated bixin or samples known to contain only bixin, pure methanol was used to redissolve the colour. If filtering was required before the sample was injected, it was necessary to ensure that the filter did not retain colour and/or reduce the colour concentration of the sample.

Results and discussion

Extractions

Hexane may be used to remove fats and oils, but petroleum ether evaporates more readily, and appears to produce a slightly cleaner extract. Drying agents, such as anhydrous sodium sulphate and Drierite®, have been used previously (Smith *et al.* 1983) to filter the chloroform extracts and prevent moisture from reaching the collection flask. However, research indicated that both these compounds can retain significant amounts of colour and adversely affect recovery. Instead, a chloroform-washed glass wool pledget was utilized.

Emulsions generated varied from sample and remained in the funnel. The emulsion was found to contain up to 17% of the total colour extracted from a cheese sample. Methanol proved to be an effective extraction solvent for emulsions. In some instances, if the methanol solution was not cloudy or somewhat viscous, it was possible to combine it directly with the chloroform extract before evaporation and HPLC, and the need for a separate LC analysis was eliminated.

Candy

The candy-making approach, utilized in this laboratory has been described previously (Lancaster and Lawrence 1986, 1992). Other researchers (Smith *et al.* 1983) described a method for the preparation of 'boiled sweets', but also reported a number of problems. For example, norbixin precipitated 'as a fine dispersion and some may have been trapped in the scum-line at the wall of the pan instead of being incorporated in the body of the batch'. By ensuring that the colour was dissolved completely in a basic aqueous solution (pH > 12), and by adding citric acid at the beginning of candy-making instead of at the end, the problem was eliminated in this laboratory. Smith *et al.* also recommended a single extraction of the candy using 50 ml 0.5% acetic acid in chloroform, but our studies, utilizing an uncoloured candy which was spiked with norbixin immediately before extraction, indicated that a single extraction removed only ca 60% of the norbixin, while a second and third produced additional recoveries of 18–30% and 6–8% respectively.

A candy to which annatto was not added was spiked with norbixin at a level of 150 μ g/g immediately before extraction. Three extractions produced a total recovery of 85.5% ($n = 2$). Therefore, 14.5% of the norbixin which survives candy-making is either not extracted, or is lost during filtration or by experimental error. Recovery of


Figure 3. HPLC chromatogram of α - and β -bixin in (a) chloroform (not quantitative) and (b) methanol (456 ng bixin per 20 μ l injection). Mobile phase, methanol-2.0% acetic acid (9:1). Detector wavelength, 500 nm. Detector range a, 0-4 a.u.s. Detector range b, 1-0 a.u.s.

norbixin from three coloured candies containing 125 μ g/g ranged from 74.2% to 89.7%. This would suggest that most of the norbixin added to the hot candy survives the manufacturing process.

It was suggested previously (Smith *et al.* 1983) that chloroform or chloroform containing 0.5% acetic acid were suitable solvents for dissolving annatto extracts immediately prior to HPLC analysis. However, our tests indicated that neither were useful for quantitative LC evaluation. Resolution was poor and peak areas were not consistent. When the chloroform was evaporated and norbixin was redissolved in methanol-2% acetic acid (9:1), resolution improved dramatically. The chromatograms also indicated that α -norbixin is more susceptible to degradation than is β -norbixin, and that geometric isomerization of the less stable α -form also occurs. Solvents were also tested with bixin, and pure methanol was found to be superior to chloroform when analysing extracts. Figure 3 shows a chromatogram of bixin, dissolved in chloroform then in methanol. The improvement in resolution and peak shape was immediately obvious.

During these studies, stock solutions of 'pure' (90%) bixin crystals were utilized. It was determined that 0.08% sodium hydroxide in methanol was a suitable solvent for dissolving the crystals and was superior to methanol. Stock solutions were diluted significantly with methanol before HPLC analysis.

Using the procedure described above, a hard candy containing norbixin (139.9 μ g/g) was prepared and extracted. Recoveries from three consecutive extractions were 53.2%, 29.4% and 7.1%, producing a total recovery of 89.7% ($n = 2$). The α - to β - ratios (7.4/1, 5.5/1 and 3.8/1 respectively) decreased as the number of extractions increased, and this suggested that α -norbixin was extracted more readily than the β -form. The ratio of α - to β -norbixin contained in the first extract closely matched the ratio (7.7/1) in the concentrated norbixin used to make the candy. Those ratios and the total recovery indicated that only minor decomposition took place when the candies were prepared.

Table 1. Annatto in commercial samples of bulk cheese, cheese spread and process cheese.

Sample	Total norbixin ^a ($\mu\text{g/g}$)	Total bixin ^a ($\mu\text{g/g}$)
Process cheese spread	1.1 ± 0.1	5.9 ± 0.2
Process cheese slices	2.5 ± 0.2	nd
Skim milk cheese	15.6 ± 0.4	nd
Canadian wine cheese	16.8 ± 0.5	5.1 ± 0.2
Cheddar, medium	13.2 ± 0.5	nd
Canadian Colby	21.2 ± 0.6	nd
Cheddar, sharp	63.8 ± 1.8	nd

nd = none detected.

^a Average of two or three determinations; three extractions/determination.

Cheese

Recovery of norbixin averaged 92.6% over a range of 1 to 110 $\mu\text{g/g}$. Evaluation of bixin-spiked samples produced recoveries which were similar to those for norbixin. For example, when bixin (22.6 $\mu\text{g/g}$) was added to an uncoloured Cheddar cheese, recovery averaged 87.6% ($n = 3$). The LC detection limit of the method was found to be ca 5 ng/g for the compounds studied.

Several commercial samples of bulk or process cheese were analysed and the results are shown in Table 1. Norbixin was present in all samples; the average level was 20.6 $\mu\text{g/g}$. In addition, bixin was found in both a cheese spread and in a wine Cheddar cheese. Although research focused on the determination of norbixin, tests indicated that any bixin present in the cheese was also extracted. Because bixin is not soluble in petroleum ether, but is soluble in ethanol-aqueous ammonia and chloroform, it is extracted at the same time and in the same manner as norbixin.

Emulsions varied from sample to sample, but the methanol wash eliminated any loss of colour. A sharp Cheddar cheese, which contained more than three times the average amount of norbixin, produced practically no emulsion, nor did the process cheese spread, which contained the lowest level of annatto. Therefore, no correlation was observed between the level of annatto in the cheese and the amount of emulsion.

Extraction with petroleum ether removed fats plus a considerable amount of yellow colour, possibly β -carotene, from both the cheese spread and the process cheese, but tests demonstrated that this colour was not annatto or a synthetic food colour, and that bixin and norbixin remained in the ethanol-aqueous ammonia layer. During HPLC analysis, bixin and norbixin retention times were confirmed using authentic internal standards.

Figure 4 shows a chromatogram of a Canadian wine cheese containing both norbixin and bixin. The chromatogram of a process cheese spread shown in figure 5 reveals that this product contains a similar amount of bixin, but that the level of norbixin is much lower than in the wine cheese.

Butter and margarine

Since the last century, annatto has been a popular colouring agent for butter. Colours such as β -carotene or annatto are usually not added during spring and summer. Our study of butter and margarine took place during September and October. Most of


Figure 4. HPLC chromatogram of a Canadian wine cheese extract dissolved in methanol. Mobile phase, methanol-2.0% acetic acid (9:1). Detector wavelength, 500 nm. Detector range, 0-1 a.u.


Figure 5. HPLC chromatogram of a process cheese spread dissolved in methanol. Mobile phase, methanol-2.0% acetic acid (9:1). Detector wavelength, 500 nm. Detector range, 0-04 a.u.

the commercial butter samples which we examined had been produced in August and early September, and did not contain significant quantities of either bixin or norbixin. In addition, a major producer reported that they no longer add colour to their products. However, successful recovery studies were conducted, and it was clearly demonstrated that the method described above was suitable for monitoring or surveillance purposes, and quantitatively determining annatto, if present in either butter or margarine.

Both uncoloured butter and margarine were spiked with known levels of bixin and norbixin. Recovery of bixin averaged $93.2 \pm 6.0\%$ ($n = 6$) over a range of 0.1 to 445 $\mu\text{g/g}$. Recovery decreased significantly above 450 $\mu\text{g/g}$ because of colour saturation and the composition of the solvents. Recovery for norbixin was almost identical to that for bixin. Geometric isomerization was again observed.

Colour manufacturers reported that annatto is used to colour margarine in addition to butter. Several commercial samples of margarine bearing various brand names were examined and none were found to contain bixin or norbixin. It was also observed that six different brands of margarine, purchased at various locations, were all prepared by the same manufacturer. All samples studied did contain 'vegetable colour', which was found to be β -carotene. Therefore, it was demonstrated that, at least in this region of Canada, annatto is not a popular colouring agent for margarine.

Table 2. Annatto in commercial samples of butter.

Sample	Total bixin ¹ (ng/g)	Total norbixin ² (ng/g)
Brand A	199 \pm 12	914 \pm 25
Brand B	12 \pm 5	nd
Brand C1	14 \pm 5	19 \pm 5
Brand C2	19 \pm 5	27 \pm 5
Brand D	17 \pm 5	33 \pm 5

nd = none detected.

¹ Average of two or three determinations;
three extractions/determination.


Figure 6. HPLC chromatogram of a commercial butter extract (containing bixin (199 ng/g) and norbixin (914 ng/g)) dissolved in methanol. Mobile phase, methanol-2.0% acetic acid (9:1). Detector wavelength, 500 nm. Detector range, 0-1 a.u's.

Five commercial samples of butter were analysed, and the results are shown in table 2. Levels of annatto in butter (ng/g range) are traditionally much lower than in cheese ($\mu\text{g/g}$ range). It is interesting to note that although it is frequently reported that bixin is the annatto compound which is added to butter, norbixin was also found in the samples that we analysed.

β -Carotene may have been added to some of the samples. However, in addition to fats and oils, β -carotene is extracted by petroleum ether and does not appear in the HPLC solution. A quantity of β -carotene was dissolved in methanol and examined using the same LC mobile phase and wavelength described above. β -Carotene was not detected at 500 nm. When the wavelength was shifted to 470 nm, a very broad peak appeared at ca 3 min. Therefore, even if β -carotene was added to butter or margarine and a trace survived the ether extraction, it did not interfere with the determination of annatto.

Some of the butter samples were found to contain trace amounts of annatto (table 2). The levels were so low in fact that one might first assume that no manufacturer would intentionally add those amounts because the colour would not be visible. One might then conclude that the peaks which appeared on the LC chromatogram were not annatto, but some other compound. However, by utilizing authentic standards, it was demonstrated that low levels of annatto were indeed present. It was determined that the source of the annatto in these products was whey cream, a by-product of the cheese-manufacturing process, which is used to make butter. Figure 6 shows a chromatogram of a commercial sample of butter, which was found to contain bixin (199 ng/g) and norbixin (914 ng/g).

References

- ANAKAWA, A., HIRATA, K., OGAWARA, T., and ÔNISHI, K., 1984, Determination of oil soluble natural dyes in foods by high-performance liquid chromatography. *Bunseki Kagaku*, 33, 586-590.
- KANADA, H., WARABI, Y., SATÔ, E., YAMASHITA, T., and TAKESHITA, R., 1981, Detection of some natural dyes by polyamide thin-layer chromatography. *Eisei Kagaku*, 27, 50-55.
- LANCASTER, F. E., and LAWRENCE, J. F., 1986, Thermal decomposition of the food colours amaranth, sunset yellow FCF and tartrazine in the presence of sucrose and glucose. *Food Additives and Contaminants*, 3, 295-304.
- LANCASTER, F. E., and LAWRENCE, J. F., 1992, Determination of total non-sulphonated aromatic amines in soft drinks and hard candies by reduction and derivatization followed by high-performance liquid chromatography. *Food Additives and Contaminants*, 9, 171-182.
- MCKEOWN, G. G., 1961, Paper chromatography of bixin and related compounds. *Journal of the Association of Official Agricultural Chemists*, 44, 347-351.
- MCKEOWN, G. G., 1963, Composition of oil-soluble annatto food colors. II. Thermal degradation of bixin. *Journal of the Association of Official Agricultural Chemists*, 46, 790-796.
- MCKEOWN, G. G., 1965, Composition of oil-soluble annatto food colors. III. Structure of the yellow pigment formed by the thermal degradation of bixin. *Journal of the Association of Official Agricultural Chemists*, 48, 835-837.
- MCKEOWN, G. G., and MARX, E., 1962, The composition of oil-soluble annatto food colors. *Journal of the Association of Official Agricultural Chemists*, 45, 761-766.
- MIYAMOTO, F., SAEKI, M., KAMUO, M., KANADA, H., NAKAGAKI, T., NISHIMURA, M., ITÔ, Y., and TAKESHITA, R., 1991, Systematic separation of artificial and natural dyes in foods and their qualitative determination by thin-layer chromatography. *Eisei Kagaku*, 37, 542-551.
- PRESTON, H. D., and RICARD, M. D., 1980, Extraction and chemistry of annatto. *Food Chemistry*, 5, 47-55.
- REITH, J. F., and GIELEN, J. W., 1971, Properties of bixin and norbixin and the composition of annatto extracts. *Journal of Food Science*, 36, 861-864.
- SMITH, P. R., BLAKE, C. J., and PORTER, D. C., 1983, Determination of added natural colours in foods. 10. Annatto. *Leatherhead Food RA, Research Report No. 431*, October (Leatherhead, UK: Leatherhead Food RA).
- TIRUMANJIA, A. S. L., 1981, Study of the carotenoid pigments of *Bixa orellana* L. seeds by TLC. *Mikrochimica Acta*, 2, 11-16.

**ANALYTICAL METHODS FOR NATURAL FOOD COLOUR
SPECIFICATIONS**

Curcumin (E100)	13
Tests	13
Colour reaction test C	13
Colour reaction test D	13
TLC test E	14
Method of assay	16
Spectrophotometry	16
HPLC	19
Annatto (E160b)	45
Definition	45
Chemical names	45
Characteristics	46
Characteristic identification tests	46
Aromatic thermal degradation products	49
Tests	50
Identification tests	50
Method of assay	50
Spectrophotometry	50
HPLC	51

2. Curcumin (E100)

The JECFA specification referred to here is that laid down in Food and Nutrition Paper 52 (Add.3) for curcumin [JECFA 1995]. Two samples of reagent grade curcumin and one sample of water-dispersible turmeric extract were used as test materials (Table 2.1.).

2.1. Tests

2.1.1. Colour reaction test C

The following tests are prescribed: "A solution of the sample in ethanol is characterised by a pure yellow colour and light green fluorescence; if this ethanol extract is added to concentrated sulfuric acid, a deep crimson colour is produced" It is recommended that this should read :

"A dilute (*ca.* 0.01% w/v) solution of curcumin in ethanol is characterised by a pure yellow colour and light green fluorescence under ultra-violet light and by its reaction with concentrated sulfuric acid; Place *ca.* 1 ml of the sample solution in a 10 ml test tube (or similar container). Using a Pasteur pipette, carefully add *ca.* 0.5-1.0 ml of concentrated sulfuric acid down the side of the inclined tube so that it forms a discrete layer at the bottom. If curcumin is present, an intense purple layer will be visible at the interface of the layers. Upon shaking (CARE), the mixture assumes a crimson colour."

2.1.2. Colour reaction test D

The following test is prescribed: "Treat an aqueous or dilute ethanolic solution of the sample with hydrochloric acid until a slightly orange colour begins to appear. Divide the mixture into 2 parts and add some boric acid powder or crystals to one portion. Marked reddening will be quickly apparent, best seen by comparison with portion to which the boric acid has not been added. The test may also be made by dipping pieces of filter paper in ethanolic solution of the colouring matter, drying at 100°C, and then

moistening with a weak solution of boric acid to which a few drops of hydrochloric acid have been added. On drying, a cherry-red colour will develop."

The following regime is recommended:

"Treat *ca.* 5 ml of a dilute ethanolic or aqueous solution of sample with *ca.* 2 ml of concentrated hydrochloric acid whereupon a permanent slight orange colour is attained. Divide the mixture into two parts and add boric acid (*ca.* 50 mg) to one (test) portion. A marked reddening of the test portion will be apparent compared to the untreated portion. The test may also be made by dipping pieces of filter paper in dilute ethanolic sample solution, drying at *ca.* 100°C and moistening with a solution of boric and oxalic acid solution* On drying, a cherry-red colour will develop."

[*0.25M boric acid (1.5 g/100 ml) + 0.25M oxalic acid (2.3 g/100 ml) in 2M hydrochloric acid.]

2.1.3. TLC test E

The following test is prescribed: "Spot 5 µl of test solution (0.01 g of sample in 1 ml of 95% ethanol) on a TLC (Microcrystalline cellulose, 0.1 mm) plate in a developing chamber containing the mixture 3-methyl-1-butanol/ethanol/water/ammonia (4:4:2:1) as solvent and allow the solvent front to ascend 10-15 cm. Examine under daylight and under UV light and observe:

- Three yellow spots with R_f between 0.2 and 0.4 under daylight and UV light
- Spots with R_f about 0.6 and 0.8 under UV light

All spots show distinct fluorescence under UV light."

Inconsistent results were found with the prescribed system using cellulose plates and 3-methylbutanol/ethanol/water/ammonia (4:4:2:1) solvent. It was not possible to distinguish three yellow spots with R_f between 0.2 and 0.4 under daylight and UV light, and spots with R_f about 0.6 and 0.8 under UV light as prescribed. Whilst three components were visible, the spots tailed badly and overlapped to such an extent as to make R_f calculations very difficult. Moreover, the main spot assumed a deep crimson/purple colour whilst developing due to the presence of ammonia which had to be removed by air-drying to reveal the prescribed yellow colour. It is conceivable that, in the presence of light, the alkaline medium was causing degradation of the curcuminoids *in situ*.

Curcuminoids can be separated on a TLC system comprising silica plates and chloroform/acetic acid (80:20) solvent [Janßen and Gole, 1984]. However, to avoid the use of chloroform, a system was developed at CSL employing dichloromethane/acetic acid where a volume mixture of 90:10 was found to give good results on silica gel-60 x 0.250 mm plates. Three discrete spots were separated with consistent R_f values in the ranges (I) 0.48-0.51, (II) 0.62-0.65 and (III) 0.77-0.81, attributable to (I) bisdesmethoxycurcumin, (II) desmethoxycurcumin and (III) curcumin respectively. Reproducible R_f values could only be obtained under carefully controlled conditions *i.e.* plate integrity, solvent mixture and full equilibration of developing chamber prior to analysis. This system was found to be more rapid than the prescribed procedure and good separation could be achieved within about 30 minutes.

When viewed under UV light the following colours were observed: (I) pale green fluorescence. (II) and (III) pale green /yellow fluorescence. Further confirmation was achieved by spraying with a mixture of 0.25M boric acid

(1.5 g/100 ml) + 0.25M oxalic acid (2.3 g/100 ml) in 2M hydrochloric acid, which was developed from a reagent described previously [Janßen and Gole, 1984]. Following treatment, the plates were heated at *ca.* 100°C for *ca.* minutes and the following respective colour changes observed :

- (I) deep red in daylight, red fluorescence under UV light
- (II) brick red in daylight, dark orange/red -fluorescence under UV light
- (III) orange/red in daylight, orange fluorescence under UV light.

These colours were stable for several hours. When this system was tested on all three curcumin / turmeric colour formulations, several minor spots were observed eluting between the main spots with signs of tailing evident in the latter. These are thought to be due to minor co-extractives of turmeric and keto-enol tautomerism respectively. This improved TLC method is therefore recommended for inclusion in the specification for curcumin.

2.2. Method of assay

2.2.1. Spectrophotometry

The JECFA specification for curcumin prescribes a spectrophotometric method of assay based on absorbance in ethanol at 425 nm and $E_{1\%}^{1\text{cm}} = 1607$. The value for the extinction coefficient is scientifically valid and generally accepted as an 'average' value for the three component curcuminoids collectively. However, several values for the wavelength of maximum absorption (in ethanol) may be found in the literature. The EC specification quotes 426 nm, whereas other sources quote between 424 and 430 nm [Hendry and Houghton, 1992; Govindajaran, 1982]. This disparity arises because of the contribution of each curcuminoid to the mixture, since each

has a different λ_{\max} value *i.e.* reported in 95% ethanol to be 430 nm (curcumin), 423 nm (desmethoxycurcumin) and 418 nm (bisdesmethoxycurcumin) [Govindajaran, 1980]. The relative proportions of each component will determine the overall λ_{\max} value in a mixture, but a value of between 424-426 nm (or expressed as 425 ± 1 nm) would seem logical since the differences in absolute absorption between these wavelengths is negligible.

A recommended sample weight of 0.08 g and a dilution factor of 20,000 are specified but if 100% dye purity is assumed, a theoretical absorbance value of 0.643 would be observed. This value is outside of the range specified in the US FDA method of 0.2 - 0.4 AU [Marmion, 1984]. This range is recommended to both minimise those errors associated with low absorbance readings and those due to departure from Beer-Lambert linearity at high absorbance values (see below). The EC specification prescribes an assay of not less than 90% total colouring matters content [EC 1995]. This would equate similarly to a theoretical absorbance value of 0.579, whereas a sample weight of 0.05 g similarly diluted would give an absorbance value of 0.362. For convenience, it would be easier to allow for an adjustment to either the final dilution or aliquot volume in order to bring the sample absorbance within the recommended range. Nevertheless, a sample weight tolerance of ± 0.001 g equates to an error in the purity value of $\pm 2.5\%$ under these conditions. When considered alongside other sources of error this is most likely to be the greatest source and should therefore be minimised by quoting the sample weight tolerance as such.

The use of modern spectrophotometers may allow the linear absorbance range to be extended beyond 0.4 AU. At CSL, a double-beam, dual source scanning spectrophotometer (Model Lambda 3, Perkin Elmer, Beaconsfield, UK - some 12 years old) was used to calibrate the response of curcumin

solutions in ethanol at 425 nm using a 1 cm quartz cell. The absorbance response was linear up to 1.7 AU ($r^2=0.9996$). Lower specification instruments would require similar calibration to determine the linear dynamic range as an alternative to specifying an absorbance range. Linear regression analysis was performed on the results and are given in Table 2.2.

For routine curcumin assays, it should be sufficient to measure the absorbance of the sample diluted appropriately in ethanol, which should be of suitable analytical grade and should be neutral. The solvent should however be checked for alkalinity or excess acidity since these may decrease colour stability rapidly. A pH of between 4 and 7 has been recommended [Govindajaran, 1980]. (Note: The concept of pH is truly meaningful only for dilute aqueous solutions. In water-like solvents such as methanol and ethanol, similar concepts may be developed [Cotton and Wilkinson, 1972]. For the purpose of this study, the pH concept should be acceptable).

Curcumin exhibits strong fluorescence in organic solvents but whilst methods for the direct measurement of curcumin have been reported the optimisation of analytical parameters have not been established. A fluorimetric method for the determination of curcumin in yoghurt and mustard using direct and derivative techniques was therefore developed [Navas Diaz and Ramos Peinado, 1992]. This method, whilst demonstrated to be simple, rapid and reliable may not be suitable for recommendation in the analysis of curcumin colouring materials since the instrumentation required may not be readily available or may prove to be too expensive for general use.

It is therefore recommended that:

- the wavelength for measurement should be standardized at either 425 or 426 nm in line with either the EC or JECFA specification

- a sample weight tolerance of 0.05 ± 0.001 g should be specified if the EC minimum colour content of 90% applies; and
- an absorbance value range of 0.2-0.4 AU should be specified and there should be appropriate adjustment of second dilution volume to bring the absorbance reading within the recommended range *unless* a wider instrumental linear dynamic range can be demonstrated.

2.2.2. HPLC

The separation and individual quantitative analysis of the three curcuminoids is important because their extinction values and stabilities vary considerably. Several methods for the HPLC analysis of curcumin have been reported which are generally straightforward. Baseline separation of the three curcuminoids was reported using an amino ($-\text{NH}_2$) bonded-phase column and ethanol mobile phase [Tønnesen and Karlsen, 1983], but the peaks were somewhat broad and tailed significantly. Nevertheless, this method was reported to give a 10-fold increase in sensitivity when fluorescence detection was used *i.e.* $\text{Ex} = 420$ nm, $\text{Em} = 470$ nm, compared to UV-VIS detection. A method which could partially separate the three curcuminoids and in turn could be separated from annatto colour components simultaneously using isocratic and non-linear gradient elution has been reported [Rouseff, 1988]. A generally-available reverse-phase system was used comprising an ODS (C_{18}) column with a tetrahydrofuran (THF):water mobile phase. THF is known to form reactive peroxides upon storage and HPLC grade THF is usually not stabilised against their formation. The peroxides may be damaging to the HPLC column and may cause sample degradation hence its use is not recommended.

An isocratic reverse-phase HPLC method developed at CSL for the analysis of annatto colouring components [Scotter *et al.*, 1994] has been successfully used for the simultaneous analysis of curcuminoids. This system comprised an HRPB C₈/C₁₈ mixed phase column (HiChrom, Reading) and an acetonitrile:aqueous acetic acid mobile phase (65:35 v/v). Baseline separation of the three curcuminoids was achieved within 5 minutes and separation of annatto components, notably 9-*cis*-norbixin which often occur together in commercial colour formulations, within 30 minute (Figure 2.1.). The curcuminoids were readily tentatively identified from their UV-VIS spectra using photodiode array (PDA) detection (Figure 2.2.). The detector response at 425 nm for each curcuminoid was linear between 0 and 45 mg/l nominal concentration. The chromatographic and spectrophotometric data obtained for the curcuminoids using this method are given in Table 2.2.

The three curcuminoids are not readily available as pure reference materials hence they either have to be prepared synthetically or purified by preparative chromatographic methods. Moreover, the three curcuminoids each have significantly different extinction values because of the relative degrees of methoxylation, hence more information on these aspects is required. The quantitative analysis of curcumin may be achieved by summing the total curcuminoids and expressing them as curcumin *i.e.* as % total colour content. This method is therefore recommended for inclusion in the specification for the assay of curcumin.

5. Annatto (E160b)

The JECFA specification originally referred to here is that laid down in Food and Nutrition Paper 52 for annatto extracts [JECFA 1992]. Since this work was completed, separate revised specifications for (i) oil- and alkali-

extracted and (ii) solvent-extracted annatto have been published in the form of a fourth Addendum [JECFA 1996]. Where the 1992 and 1996 specifications differ with respect to those criteria which have been evaluated here, appropriate comments have been made. Eight commercially-available annatto formulations comprising both bixin and norbixin types were obtained for testing and analysis (Table 5.1.). Where appropriate, purified samples of bixin and norbixin were prepared in the laboratory.

5.1. Definition

5.1.1. Chemical names

Inconsistencies in the systematic naming of annatto colour principles in the specification have been identified [Annex 1]. Rules for the nomenclature of carotenoids using a semi-systematic approach were approved by IUPAC in 1974 and have been reproduced recently [Weedon and Moss, 1995]. The carotenoid stem name implies *trans*- configuration unless the contrary is indicated. Following the designation of absolute configuration (if any), geometrical configuration is indicated by citing the double bond or bonds with a *cis*- configuration. The following nomenclature is recommended :

trans-bixin : methyl hydrogen 6,6'-diapocarotene-6,6'-dioate

'natural' (*cis*-) bixin : 6'-methyl hydrogen 9'-*cis*-6,6'-diapocarotene-6,6'-dioate

trans-norbixin : 6,6'-diapocarotene-6,6'-dioic acid

'natural' (*cis*-) norbixin : 9'-*cis*-6,6'-diapocarotene-6,6'-dioic acid

5.2. Characteristics

5.2.1. Characteristic identification tests

The characteristic identification Test (B) Spectrophotometry, has been revised in the 1996 specification. The test has been redefined as "Test B. UV/VIS absorption." Different criteria have been applied to oil- and alkali-extracted annatto, and to solvent-extracted annatto preparations. In the specification for oil- and alkali-extracted annatto the test states that:

"Oil-extracted annatto, diluted with acetone shows an absorbance maximum at about 439, 470 and 501 nm." and

"Alkali-extracted annatto, diluted with water, shows an absorbance maximum at 451-455 nm and 480-484 nm."

Thus acetone has replaced chloroform as the solvent, the wavelength criteria for oil-extracted annatto are approximate and the wavelength criteria for alkali-extracted annatto are now given as ranges. For solvent-extracted bixin preparations the test criteria differ *i.e.* the following are stated:

Bixin preparations in acetone show an absorbance maximum at about 502 nm." and -

"Norbixin preparations in 0.1N potassium hydroxide solution show an absorbance maximum at about 482 nm."

Similarly, acetone has replaced chloroform as the solvent for bixin preparations and 0.1N potassium hydroxide is specified for norbixin preparations in place of water. The wavelength criteria for both types of preparation are approximate.

The reasons for the differences in the test criteria between oil- and alkali-extracted annatto, and solvent-extracted preparations are not clear. The overall shape and structure of carotenoid spectra is diagnostic and generally effects the degree of planarity that the chromophore can achieve. The degree of fine structure decreases when coplanarity is lost due to steric strain which can be due, amongst other things, to the presence of conjugated carbonyl groups, as is the case for bixin and norbixin. Hence they exhibit typical carotenoid UV-VIS spectra characterised by two distinct central absorption maxima separated by *ca.* 30 nm and two less distinct peripheral maxima, usually present as shoulders or points of inflexion. The relative intensity and position (and hence specific absorbance) of the main maxima are determined largely by geometric configuration and the solvent used and have been discussed in detail elsewhere [Scotter *et al.*, 1994].

Whilst 9'-*cis*-bixin is the naturally-occurring form, the various annatto extraction processes promote isomerisation to other mono-*cis* (and to a lesser extent di- and tri-*cis*) isomers in addition to the *trans* isomer. Moreover, extraction into aqueous alkali is effected by hydrolysis to the corresponding isomers of norbixin, and any thermal processes provide the potential for degradation with subsequent formation of yellow-coloured conjugated species and aromatic products [Scotter, 1995]. The combination of such processes generally results in the production of annatto formulations which comprise mixtures of isomers. These are usually characterised by the presence of one or two major isomer *i.e.* 9'-*cis*- and /or *trans*- (bixin) or 9'-*cis*- (norbixin), plus lower levels of minor *cis*-isomers. *Trans*-norbixin is sparingly soluble in most solvents and therefore does not contribute significantly to the isomer profile of water-soluble forms of annatto.

However, its presence can provide qualitative information relevant to formulation production history. The spectrophotometric characteristics of annatto formulations are therefore generally less distinct than those of the corresponding pure compounds, especially when there are relatively large amounts of non-bixin/norbixin coloured degradation products present. Norbixin is particularly susceptible to oxidation in air hence norbixin formulations often have relatively low norbixin content. This can lead to severe spectrophotometric distortion characterised by changes in relative absorption intensities rather than by noticeable shifts in λ_{\max} values. The changes made in the 1996 specifications appear to go some way to address these aspects *i.e.* by specifying approximate λ_{\max} values or wavelength ranges. Whilst the improved characteristic UV/VIS identification tests for the different types of annatto preparations are recommended, more information on the reasons for the difference in specification criteria between oil-/alkali-extracted and solvent extracted annatto preparations is required.

Test (D) TLC. Following discussions between CSL and MAFF it was agreed that TLC procedures for annatto would not be investigated. Moreover, this test has been removed from the revised specifications for annatto [JECFA 1996].

Test (E) Colour reaction. Both bixin and norbixin give a characteristic dark blue colour reaction with concentrated sulfuric acid. Six different annatto formulations were tested using the following procedure: a small amount of sample (*ca.* 250 mg) was placed in a 10 ml glass test tube and *ca.* 1 ml of concentrated H_2SO_4 added carefully down the inclined side of the tube. The tube contents were mixed by swirling and any colour change observed. The results are given in Table 5.3 and show that the test appears to be reliable only for relatively pure solid formulations or those which are essentially

free from interfering co-formulation materials. It was recommended that due to limitations in scope, this test be excluded from revised specifications for annatto. However, this test has not been included in the revised specifications for annatto [JECFA 1996].

5.2.3. Aromatic thermal degradation products

Studies on characterisation of the coloured thermal degradation products formed during the commercial extraction of the food colouring material annatto [Scotter, 1995], have shown that the major degradation product of the principal colouring agent 9'-*cis*-bixin is the *trans*- monomethyl ester of 4,8-dimethyltetradecahexaenedioic acid (C_{17}). The formation of C_{17} was shown to be accompanied by the formation of *m*-xylene. Whilst thermodynamic considerations indicated that elimination schemes involving scission of the 9'-*cis*- bond were most likely to predominate, it was postulated that other competitive elimination processes could occur. These processes could lead to the formation of analogous C_{18} and C_{13} products with the corresponding formation of toluene and dimethyl dihydronaphthalene, all of which are undesirable in colouring materials intended for food use. In order to aid in the establishment of specific purity criteria in support of regulations [EC, 1994], suitable methods of analysis were required to determine the presence and levels of aromatic degradation products in commercial annatto formulations. Such a method has been developed which uses ambient alkaline hydrolysis, followed by solvent extraction and capillary gas chromatography. Twenty samples of commercially-available annatto formulations have been analysed for toluene and *m*-xylene. Fifteen of the 20 samples analysed contained < 5 mg/kg, 4 samples contained between 5 and 10 mg/kg and one sample contained 12 mg/kg of toluene. Seven samples contained < 5 mg/kg and a further 7 samples contained between 7 and 25 mg/kg *m*-xylene. Of the remaining 6 samples, four contained between 30 and 88 mg/kg and two

contained 160 and 200 mg/kg *m*-xylene. The results provide evidence for the thermal degradation of annatto during processing resulting in the contamination of both bixin and norbixin annatto formulation types with *m*-xylene. Two samples of norbixin of known production history were analysed specifically to identify possible differences in their degradation component profiles. They were found to differ significantly in their *m*-xylene contents which is consistent with their respective production histories [Annex 2]. This method is therefore **recommended** for inclusion in the purity specification for all annatto types.

5.3. Tests

5.3.1. Identification tests

Test (C) Column chromatography and Carr-Price reaction. This test is based on the characteristic colour reaction (orange-red to blue green) between bixin (and norbixin) with Carr-Price reagent, following column chromatographic purification on alumina. The main problem identified with this test was with the use of benzene as solvent, hence it was **recommended** that a suitable alternative solvent should be found [Annex 1]. In the revised [JECFA 1996] specification for annatto however toluene is specified instead.

5.4. Method of assay

5.4.1. Spectrophotometry

The assay procedures prescribed in the 1992 JECFA specifications have been modified in the revised specifications [JECFA 1996]. For the assay, quantitative analysis is achieved by reference to the extinction coefficient ($E_{1\text{cm}}^{1\%}$) values for both bixin and norbixin, either in acetone (which has replaced chloroform) or in "0.1N" sodium hydroxide respectively. The extinction coefficient and λ_{max} values are identical to those prescribed in the

EC specification for annatto, where chloroform and "KOH solution" are the prescribed solvents for bixin and norbixin preparations respectively.

The $E_{1\text{cm}}^{1\%}$ value of 2870 given in the EC specification however appears to be an industry-derived 'average' and has no relevance to literature values [Scotter *et al.*, 1994], which are derived from pure substances (Table 5.2.). Whilst it would be preferable from an analytical viewpoint to improve the assay by inclusion of the data given, the assay method would become more complicated depending upon the formulation type analysed [Annex 3]. Therefore the recommendation is that if the assay requires a single λ_{max} and $E_{1\text{cm}}^{1\%}$ value, then the results *must* be expressed in terms of the major isomer present *i.e.* 9'-*cis*- rather than "bixin" or "norbixin". Such procedures are limited by the presence of non-bixin/-norbixin absorbing species as discussed above (5.2.1.), hence an alternative assay procedure using HPLC has been developed (5.4.2.).

5.4.2. HPLC

HPLC has been used successfully for the analysis of annatto formulations [Scotter *et al.*, 1994; Scotter, 1995] in which a simple isocratic elution system was used to separate not only the major and minor colouring principals, but also several higher *cis*-isomers and the main coloured thermal degradation product, C_{17} . In a previous study, 21 different bixin and norbixin formulations including solid materials, aqueous and oil solutions, and spray-dried / encapsulated formulations were analysed using this method and the results compared with those obtained using spectrophotometric methods [Annex 3]. HPLC was shown to provide qualitative and quantitative data, particularly with respect to the determination of coloured degradation compounds. A typical separation is given in Figure 5.1. The system comprised an HRPB C_8/C_{18} column (HiChrom, Reading) and an isocratic acetonitrile : aqueous acetic acid

mobile phase (65:35 v/v). PDA detection at 450 nm with 10 nm bandwidth (for monitoring) was used to characterise the separated components by their spectra. Separation of the major and minor colour principles and coloured degradation products was achieved within 30 minutes.

The 9'-*cis*- and *trans*- isomers of bixin and norbixin, and *trans*-C₁₇ are not commercially-available as pure reference materials and require laboratory preparation and purification [Scotter *et al.*, 1994]. It is recommended that the total amount of bixin (or norbixin) isomers should be expressed in terms of the major 9'-*cis*-isomer, hence reference materials for these isomers need to be prepared. Bixin is available commercially as a '95%' pure crystalline product from which 9'-*cis*-bixin may be readily purified by crystallisation, and 9'-*cis*-norbixin by alkaline hydrolysis followed by crystallisation. Bixin (or norbixin) content and isomeric purity can then be assessed using spectrophotometry and HPLC respectively. It should be noted that even as isolated solids these compounds are susceptible to oxidative degradation (especially norbixin), and to isomerisation once in solution. Conditions of storage are therefore important. It is therefore recommended that HPLC be included in the specifications as a suitable alternative assay method for annatto.

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